

## SINGLE WAVELENGTH, X-SERIES

### Long wave, LW (UV-A, 365nm)

- Integrally filtered tube (BLB) or unfiltered tube (BL)
- With or without filter

Model	Equipped with	Typical Peak UV Intensity ( $\mu\text{W}/\text{cm}^2$ ) at 10" (25cm)
X-15A	One 15-Watt BLB, BLE-1800B	800
X-15N	One 15-Watt BL, BLE-1T151 with filter	600
X-30	One 30-Watt BLB, BLE-5000B	1,100
X-40	One 40-Watt BLB, BLE-7900B	1,400
XX-15A	Two 15-Watt BLB, BLE-1800B	1,100
XX-15N	Two 15-Watt BL, BLE-1T151 with filter	850
XX-40	Two 40-Watt BLB, BLE-7900B	1,850
XX-40N	Two 40-Watt BLB, BLE-7900B with filter	1,100

### Medium wave, MW (UV-B, 312nm)

- Unfiltered tube
- With LONGLIFE filter

Model	Equipped with	Typical Peak UV Intensity ( $\mu\text{W}/\text{cm}^2$ ) at 10" (25cm)
X-15B	One 15-Watt MW, BLE-1T158	660
XX-15B	Two 15-Watt MW, BLE-1T158	1,150

### Short wave, SW (UV-C, 254nm)

- Unfiltered tube
- With LONGLIFE filter

Model	Equipped with	Typical Peak UV Intensity ( $\mu\text{W}/\text{cm}^2$ ) at 10" (25cm)
X-15F	One 15-Watt SW, BLE-1T155	640
XX-15F	Two 15-Watt SW, BLE-1T155	1,100

## MULTI WAVELENGTH, X-SERIES

- Unfiltered long, medium and short wave tubes
- With LONGLIFE filter

### Long wave (UV-A, 365nm) / Short wave (UV-C, 254nm)

Model	Equipped with	Typical Peak UV Intensity
-------	---------------	---------------------------

		( $\mu\text{W}/\text{cm}^2$ ) at 10" (25cm)
XX-15NF	One 15-Watt BL, BLE-1T151 One 15-watt SW, BLE-1T155	UV-A: 550 UV-C: 490

### Long wave (UV-A, 365nm) / Medium wave (UV-B, 312nm)

Model	Equipped with	Typical Peak UV Intensity ( $\mu\text{W}/\text{cm}^2$ ) at 10" (25cm)
XX-15NB	One 15-Watt BL, BLE-1T151 One 15-watt MW, BLE-1T158	UV-A: 550 UV-B: 590

### Medium wave (UV-B, 312nm) / Short wave (UV-C, 254nm)

Model	Equipped with	Typical Peak UV Intensity ( $\mu\text{W}/\text{cm}^2$ ) at 10" (25cm)
XX-15BF	One 15-Watt MW, BLE-1T158 One 15-watt SW, BLE-1T155	UV-B: 590 UV-C: 490

### GERMICIDAL, X-SERIES

- Unfiltered short wave (UV-C) tubes without filter

Model	Equipped with	Typical Peak UV Intensity ( $\mu\text{W}/\text{cm}^2$ ) at 10" (25cm)
X-15G	One 15-Watt SW, BLE-1T155	820
X-30G	One 30-Watt SW, BLE-1T305	1,850
XX-15G	Two 15-Watt SW, BLE-1T155	1,550