

B-260 LONG WAVE BATTERY LAMP - SPECIFICATION

Tubes (2)	6-watt UV-A Tube (365nm), Part No. BLE-480B
Power Supply	6-volt, 4 amp-hour Rechargeable Battery
Charger	Built-in AC Recharger (230 Volt)
Charging Time	Approximately 24 hours
Operating Time (full charge)	Dual Tubes : Approximately 5 hours
Dimensions (H x W x D)	10.25 x 3.15 x 5.63 in (26 x 8 x 14.3 cm)
Weight	4 lbs. (1.8 kg)