

Spectroline E Series UV Lamps - Specification

Product Number	UV [nm]	Type	W x L x H cm	Replace tube	Replace filter
EA-140	365	1 x 4w BLB tube	8.5 x 18 x 6.5	BLE-220B	
EN-140	365	1 x 4 watt tube	8.5 x 18 x 6.5	BLE-270W	
EN-140L	365	1 x 4w tube with filter assembly	8.5 x 20.5 x 6.5	BLE-270W	2F003
EA-240	365	2 x 4w BLB tube	8.5 x 18 x 6.5	BLE-220B	
EA-150	365	1 x 5w BLB tube	8.5 x 20.5 x 6.5	BLE-50H	
EA-160	365	1 x 6w BLB tube	8.5 x 26 x 6.5	BLE-480B	
EN-160	365	1 x 6w tube	8.5 x 26 x 6.5	BLE-6T365	
EN-160L	365	1 x 6w tube with filter assembly	8.5 x 26 x 6.5	BLE-6T365	2F005
EA-180	365	1 x 8w BLB tube	8.5 x 34 x 6.5	BLE-760B	
EN-180	365	1 x 8w tube	8.5 x 34 x 6.5	BLE-8T365	
EN-180L	365	1 x 8w tube with filter assembly	8.5 x 34 x 6.5	BLE-8T365	2F082
EN-280L	365	2 x 8w tubes with filter assembly	8.5 x 34 x 6.5	BLE-8T365	2F082
EB-160C	312	1 x 6w tube with filter assembly	8.5 x 26 x 6.5	BLE-6T312	2F006
EB-180C	312	1 x 8w tube with filter assembly	8.5 x 34 x 6.5	BLE-8T312	2F018
EB-280C	312	2 x 8w tube with filter assembly	8.5 x 34 x 6.5	BLE-8T312	2F018
EF-140C	254	1 x 4w tube with filter assembly	8.5 x 20.5 x 6.5	BLE-2537S	2F001
EF-160C	254	1 x 6w tube with filter assembly	8.5 x 26 x 6.5	BLE-6254S	2F006
EF-180C	254	1 x 8w tube with filter assembly	8.5 x 34 x 6.5	BLE-8T254	2F018
EF-260C	254	2 x 6w tube with filter assembly	8.5 x 26 x 6.5	BLE-6254S	2F006
EF-280C	254	2 x 8w tube with filter assembly	8.5 x 34 x 6.5	BLE-8T254	2F018
ENF-240C	365/254	1 x 4w 365nm, 1 x 4w 254nm with filter assembly	8.5 x 20.5 x 6.5	BLE-270W BLE-2537S	2F001
ENF-260C	365/254	1 x 6w 365nm, 1 x 6w 254nm with filter assembly	8.5 x 26 x 6.5	BLE-6T365 BLE-6254S	2F006
ENF-280C	365/254	1 x 8w 365nm, 1 x 8w 254nm with filter assembly	8.5 x 34 x 6.5	BLE-8T365 BLE-8T254	2F018
ENB-260C	365/312	1 x 6w 365nm, 1 x 6w 254nm with filter assembly	8.5 x 26 x 6.5	BLE-6T365 BLE-6T312	2F006
ENB-280C	365/312	1 x 8w 365nm, 1 x 8w 254nm with filter assembly	8.5 x 34 x 6.5	BLE-8T365 BLE-8T312	2F018
EBF-260C	312/254	1 x 6w 312nm, 1 x 6w 254nm with filter assembly	8.5 x 26 x 6.5	BLE-6T312 BLE-6T254S	2F006
EBF-280C	312/254	1 x 8w 312nm, 1 x 8w 254nm with filter assembly	8.5 x 34 x 6.5	BLE-8T312 BLE-8T254	2F018
EF-140	254	1 x 4w tube, germicidal	8.5 x 18 x 6.5	BLE-2537S	
EF-160	254	1 x 6w tube, germicidal	8.5 x 26 x 6.5	BLE-6254S	
EF-180	254	1 x 8w tube, germicidal	8.5 x 34 x 6.5	BLE-8T254	

For more information on the Spectroline E-Series UV Lamp range, visit our UVMAN website at uvman.co.uk

