

MiniMAX Series UV Lamps - Specifications				
Model	Wavelength nm	Weight g	Dimension cm	Features
UV-5A	365	290	22.9 x 5.7 x 2.5	Long wave BLB
UV-5F	254	312	22.9 x 5.7 x 3.2	Short wave tube with filter
UV-5NF	365/254	340	22.9 x 5.7 x 3.2	Long/short multiband tube with filter and sliding wavelength selector
UV-5G	254	312	22.9 x 5.7 x 3.2	Short wave tube – no filter
UV-4AW	365/white	290	20.9 x 5.7 x 3.2	Long wave BLB and white light
UV-4FW	254/white	290	20.9 x 5.7 x 3.2	Short wave tube with filter and white light
UV-4NFW	365/white	290	20.9 x 5.7 x 3.2	Long/short multiband tube with sliding wavelength selector and white light